

# Mechanical column scale with BMI display and calculation seca 755



**new**



- ◀ Quick and easy ascertainment of BMI . Simultaneous evaluation of nutritional condition . Clearly designed dial and low platform

# The intelligent all-round solution

This dial column scale can do something other mechanical scales can't: determine the Body Mass Index (BMI) at the same time as weighing. The BMI is recognised by the World Health Organisation as being the indicator of one's individual nutritional condition. When fitted with the measuring rod seca 220, the dial scale becomes an all-in-one mechanical weighing and measuring station. The elegant all-round solution is especially suitable for use in the general medical sector in doctors' surgeries and hospitals; at the same time, it is a comfortable instrument for use by healthworkers all around the world. When installed in pharmacies and fitness studios, it provides an exclusive additional service for the customers.

## Quick and easy ascertainment and evaluation of BMI


Underweight, normal weight, overweight or adiposity? For the first time, the seca 755 answers this question mechanically: When the patient mounts the scale, a "BMI disc" is set in motion just below the very precise weight display in the large dial. The disc is visible in a small window next to which a height scale is printed. The BMI can now be quite easily read off from this height scale. The coloured sectors of the BMI disc provide information on how the ascertained BMI is to be interpreted, on whether the patient is underweight, overweight or of normal weight. If there is uncertainty as regards the patient's height, this can be determined precisely in just one step with the optional integrable telescopic measuring rod seca 220. It is hardly possible to be quicker and more precise.

## Clearly designed dial and low platform

When designing this elegant scale, special attention was paid to the ergonomics. The dial which is tilted slightly towards the observer is particularly large and clearly laid out, making it easier to read. The seca 755 is easy to mount as the platform is very low. The tread area is so large that even tall and heavyweight patients have sufficient space. The tread area is also covered with a non-slip rubber mat for safe footing.

## Mobile with integrated transport castors

Thanks to the integrated transport castors, the seca 755 is easy to move around and has no problems with uneven surfaces due to its adjustable levelling feet.



## Technical data

### seca 755

Capacity: 160 kg / 350 lbs

Graduation: 500 g / 1 lb

Dimensions (WxHxD): 360 x 1,000 x 628 mm

Platform (WxHxD): 350 x 84 x 360 mm

Weight: 18 kg

Function: BMI

Optional: Measuring rod seca 220